



## Solid State Rotors

Revolution NMR supplies a wide range of rotors for solid state NMR along with drive caps and sample spacers. Our standard zirconia models include:

Outer Diameter (mm) TW = thin wall	Compatibility: Revolution (R) Chemagnetics (C) Varian® (V) Agilent® (A) Bruker®	Sample Volume (µl)	Rated Spinning Speed (kHz)
1.6*	A, V, R	8	40
2.0*	R	11	35
2.5*	B	12	35
3.2*	A, V, C, R	11	25
3.2*	A, V, C, R	22	25
3.2TW*	A, V, C, R	36	18
3.2*	B	30	24
3.2TW*	B	42	24
4.0	A, V, C, R	52	18
4.0 TW	A, V, C, R	83	15
4.0 NB	A, V, C, R	46	18
4.0	B	80	15
5.0	A, V, C, R	230	12
5.0TW	V, C, R	160	9
5.0NB	A, V, C, R	90	13
6.0	A, V, C, R	240	10
6.0NB	A, V, C, R	155	10
7.0	R	440	8
7.0	B	440	7
7.5	V, C	450	7
9.5	V, C	750	5
14TW	V, C	2100	3.0
14	V, C	2500	3.5

\*All 1.6, 2.0 and 3.2mm rotors are reduced Al

Please contact us for additional information and quotations and ask us about custom rotors

Mailing Address: 4921 Eagle Lake Drive, Fort Collins, CO 80524 P (970) 472-0613 F (970) 416-8896

Shipping Address: 2416 Donella Court, Unit B, Fort Collins CO 80524 P (970) 493-6600

[www.revolutionnmr.com](http://www.revolutionnmr.com)